

FCM Lysing Solution (Fixative Free) 10×

I.

: LYS01, LYS02, LYS03 : 100T, 250T, 500T
 : 2 - 8°C :

II.

10

III.

1. 10× 1×
- 2.
- 3.
- 4.
5. 2 ml 1× FCM Lysing Solution
6. 15
7. 300 - 400 × g 5 2 ml PBS 300 - 400 × g 5
8. 0.5 ml PBS
9. 2

IV.

- 1.
- 2.
- 3.

V.

S1001	Flow Cytometry Staining Buffer	125 ml
LSC01	FCM Lysing solution for BC (ready-to-use)	100T
LSC02	FCM Lysing solution for BC (ready-to-use)	250T
LSB01	Lysing solution for FACS 10×	100T
LSB02	Lysing solution for FACS 10×	250T
LSB03	Lysing solution for FACS 10×	500T
NGS010	Normal Goat Serum	10 ml
NRS010	Normal Rabbit Serum	10 ml
NRS500	Normal Rat Serum	500 1



MultiSciences Biotech Co., Ltd

Technical Service: tech@liankebio.com

To place an order: order@liankebio.com